

# X640 L.S.O. FITTING AND CARE INSTRUCTIONS

1. Place handles as close to the closing system as possible.
2. Be sure handles are next to closing system. Stand holding product, side panels parallel to floor. Pull outward to fully extend cord.
3. Place closing system around body, making sure closing system is centered on spine. Supports holding closing system must be equally spaced on each side of spine. Adjust up or down spine to determine best support location.
4. Once location is determined, pull slightly outward on side panels to maximize open position before using fabric fasteners to close front, right over left.
5. Pull handles away from fabric fastener and across body, placing handles anywhere on fabric fastener strip encircling side panels once desired support is achieved. Loosen the same way.

## User Tips

Pull handles only slightly at first until you get used to the power.

Tighten to stand, loosen to sit. You control the tension.

Pulling in stomach muscles before tightening makes pulling handles easier.

Do not wear tight for long periods of time unless instructed by physician.

The product was designed to enable adjustment throughout the day to provide more or less support as needed.

To disconnect closing system, simply pull closing system away from fabric fastener. Attach closing system by placing closing system on fabric fastener.

## Removal Of Anterior Panel to the LaceIT X640

1. Locate the inside lining of the left side panel.
2. Locate the anterior panel pocket.
3. Lift up the anterior panel pocket, and locate and remove the rigid plastic insert.

## Attaching the Solid Posterior Frame

1. Lift up the elastic straps located on the inside, posterior portion of the L.S.O.
2. Place the L.S.O. around the outside portion of the Solid Posterior Frame, making sure the L.S.O. superior edge is facing up.
3. Align the L.S.O. elastic straps over the Solid Posterior Frame alignment tabs, and secure the L.S.O. by sliding the elastic straps through the Solid Posterior Frame tabs.

